SUNSEPT

SUN280

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EN SPRAY AND WIPE DISINFECTANT PLUS

Ready-for-use alcoholic fast disinfectant for the disinfection of small alcohol-resistant areas.

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- ✓ QAV-free
- ✓ Aldehyde-free, phenol-free
- ✓ Tuberculocide
- ✓ Limited virucidal PLUS
- ✓ Also as impregnating fluid for wipes

DANGER
GEFAHR
DANGER
PERICOLO

 PRISMAN GmbH
Otto-Hahn-Ring 6-18
64653 Lorsch
+49 (0) 6251-866 980-0

CE 0124

BAG Nr. CHZN3442
BAuA Nr. N-77270

LOT

REF SUN280L

REF SUN280N

5 ltr. e
10 ltr. e

DE **Type of application:** As per the Medical Devices Law (MPG): Ready-for-use fast alcoholic disinfectant for the prophylactic surface disinfection of alcohol-resistant medical products for medical purposes. Not suitable for the terminal disinfection of invasive medical products. As per the Biocidal Products Directive: Ready-for-use alcoholic wet wipes for the disinfection of small alcohol-resistant areas. For disinfection, moisten a suitable cloth with the disinfectant solution and wipe the area, or apply the solution directly to the surface and then wipe. For spray disinfection, the product is sprayed on the surface to be disinfected and wiped. Make sure that the surface to be disinfected is fully covered. Then let it take effect, in line with the application times. Only for use on alcohol-tolerant materials (e.g. not suitable for acrylic glass). The disinfectant solution is only allowed to be used to soak cloths of PRISMAN GmbH, which are compatible with the solution, or with cloth materials and application times which are compatible as per EN 16615. **Application times:** Bacteria and candida albicans (VAH, no mechanical action): Conc. 30 sec. (clear conditions) / Conc. 1 min. (dirty conditions). Bacteria and candida albicans as per EN 13277, EN 13624 (dirty conditions): Conc. 30 sec. Bacteria and candida albicans as per EN 1276, EN 1650 (dirty conditions): Conc. 1 min. TdB (tuberculosis) inactivation (as per DGHM), EN 14348: Conc. 30 sec. (clean conditions). Conc. 1 min. (dirty conditions). EN 14348). Enveloped viruses (e.g. HBV, HIV and HCV) as per RKI/DVV with continuation: Conc. 30 sec. Adenoviruses (EN 14476): Conc. 30 sec. (dirty conditions). Norovirus (EN 14476): Conc. 1 min. (clean conditions). Rota viruses (EN 14476): Conc. 30 sec. (dirty conditions). "Limited virucidal effect PLUS": Conc. 1 min. Vaccinia viruses, adenoviruses and noroviruses in a carrier test as per DVV 2012: Conc. 2 min. (clean conditions). **Application times for cloth-wipe systems:** Bacteria and candida albicans as per EN 16615, VAH (with mechanical action), dirty conditions: Conc. 1 min. Bacteria and candida albicans: Entry in the VAH list: Conc. 5 minutes (shortest application time for surface disinfectants in the VAH list). **Directive 648/2004/EC:** Contains scent. **100g contains:** 37,0g 1-propanol, 24,0g ethanol. **Safety guidelines:** H 318 Causes severe eye damage. H 226 Flammable liquid and vapour. H 336 Can cause drowsiness or dizziness. Keep away from heat/sparks/naked flames/heat surfaces. No smoking. Wear eye/face protection. **IN THE EVENT OF CONTACT WITH THE EYES:** Rinse carefully for a long time with water. Remove any contact lenses if possible. Continue to rinse. If eye irritation persists: Get medical advice/attention. Seal the containers and store them in a well-ventilated place. Dispose of the contents/container in accordance with local/regional/national/international regulations. Flammable! Do not spray into naked faces! Flush spot as per DIN 51755: approx. 25 °C. Application max. 50ml per m². The total applied amount per room must not exceed 100ml per m² of room area. Keep the containers tightly sealed. Only put the packaging into the recycling system if it is empty of residue. Dispose of residual product in accordance with EAK 07 06 99.

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